



ROTHBAND

PATIENT POSITIONING

INSTRUCTIONS FOR USE

IMAGING STEPS – S103, S105

VERSION 2.0

**UK
CA**

Table of Contents

INTRODUCTION	3
INTENDED USE	3
COMPLIANCE	3
TECHNICAL INFORMATION	4
ASSEMBLY INSTRUCTIONS	4
Attaching the Handrail.....	4
CAUTIONS AND WARNINGS	5
CLEANING INSTRUCTIONS	6
CONTACT DETAILS	6

INTRODUCTION

These Instructions cover the following models:

PRODUCT
S103 - Single Step
S105 - Double Step
S107 - Handrail

This product is available as a standalone product, with a single handrail or double handrail.

INTENDED USE

The Rothband Imaging Step range has been designed for use during x-ray examinations of the lower limb.

The Rothband Single Step allows the patient to be positioned for weight bearing radiographic examinations.

As all Radiographers are aware, there are many variations in technique both departmentally and due to clinical situation and patient build. These instructions are intended as a guide to enable a local protocol to be established and may be refined and adapted as and when required.

COMPLIANCE

The Standing foot plate range complies with the Medical Device Directive (MDD) 93/42/EEC. and has been designed and manufactured to comply with the following standards:

Standard	Description
BS EN 14183:2003	Step Stools
EN ISO 14971:2012	Medical devices - Application of risk management to medical devices (ISO 14971:2007, Corrected version 2007-10-01)
EN ISO 15223-1:2016	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements.

TECHNICAL INFORMATION


PRODUCT	DIMENSIONS (LxWxH)	WEIGHT	STEP HEIGHT	WEIGHT LIMIT
S103 - Single Step	500 mm x 400 mm x 1050 mm	4 Kg	210 mm	150 Kg
S105 - Double Step	510 mm x 410 mm x 1050 mm	5 Kg	Step 1 – 210mm Step 1 - 2 – 200mm	150 Kg
S107 - Handrail	400mm x 1050 mm	1.5 Kg		N/A

ASSEMBLY INSTRUCTIONS

Attaching the Handrail

For both S103 and S105

Mount each side handle using the bolts and shakeproof washers provided. There are two bolts per handle. An Allen key is provided for the bolts.



Fix the cross piece using 4 Allen bolts and shakeproof washers.



CAUTIONS AND WARNINGS

- Individuals embarking the step should be closely supervised and assisted.
- Individuals should not be left un-supervised whilst on the step.
- There are no special skills required to use this positioning aid.
- The step is only suitable for indoor use.
- The step should only be used on a flat and level surface.
- The medical professional is responsible for applying his/her best medical judgment when positioning a patient
- The step is suitable for continuous use.
- Not suitable for sterilisation.
- Do not place device on or near a heat source.
- Intended for professional healthcare facility environments where operators with medical training are continually available when patients are present.
- Do not use abrasive cleaners, please see relevant section for the correct cleaning instructions.
- The step should be regularly inspected for damage and loose fasteners.
- The above warnings, cautions and any safety considerations should be observed on a routine and regular basis, not only upon installation.
- The device has a weight limit of 150Kg

CLEANING INSTRUCTIONS

Use a solution of mild soap or detergent or disinfectant* in hand-hot water.

Apply using a soft cloth & light pressure to prevent swirls.

Rinse off with clean water and dry by blotting with a soft cloth/towel.

WSR Medical Solutions Limited recommend Envirolex hard surface cleaner.

*DO NOT Use Solvents or Alcohol based cleaning products.

DO NOT Use any scouring materials/pads to clean the product as this will damage the surface finish.



Cross contamination is a major hazard to patients. The Surgical Patient Supports should therefore be cleaned in between patients.

CONTACT DETAILS



WSR Medical Solutions
(Trading as Rothband)
Unit 2 Hargreaves Mill
Hargreaves Street
Haslingden
Lancashire
BB4 5RQ
United Kingdom



Eurolink
Europe
Compliance
Limited,
Dublin
D02 AY86

Tel: 01706 830086

e-mail: sales@rothband.com