


Measuring for a Lead Apron - Male

Copy and Send

Name
Hospital
Phone
E-Mail


Male Sizing Guide					
Size	S	M	L	XL	XXL
Chest (cm)	82-90	92-105	110-115	120-125	130-135
Chest (Inches)	32 - 35	36 - 41	42 - 45	46 - 48	51-53
Waist (cm)	up to 78	82-88	90-100	102-113	115-125
Waist (inches)	up to 31	32-34	35-39	40-44	45-49
Apron Length					
One Piece	80 cm	85 cm	90 cm	100 cm	110 cm
Two Piece					
Top	58 cm	58 cm	58 cm	58 cm	65 cm
Kilt	53 cm	53 cm	53 cm	53 cm	55 cm
Up to 5cm addition or subtraction in length at no extra cost.					



Measuring for a Lead Apron - Female

Copy and Send

Name
Hospital
Phone
E-Mail



Female Sizing Guide					
Size	Petite	S	M	L	XL
Chest (cm)	81-86	87-96	97-106	107-118	119-130
Chest (Inches)	32-34	34-38	38-42	42-46	47-51
Waist (cm)	64-69	70-79	80-89	90-101	102-113
Waist (inches)	25-27	28-31	31-35	35-40	40-44
Hips (cm)	89-94	95-104	105-114	115-126	127-138
Hips (inches)	35-37	37-41	41-45	45-50	50-54
UK size	8-10	12-14	16-18	20-22	24-26
Apron Length					
One Piece	80 cm	85 cm	90 cm	100 cm	110 cm
Two Piece					
Top	55 cm	58 cm	58 cm	58 cm	58 cm
Kilt	50 cm	53 cm	53 cm	53 cm	53 cm